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```
// atomic number and the number of each species
```

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```
// {1:Flat Plate (2  $\theta$  only) 2:Cylindrical 3:Flat Plate ( $\theta$  -2  $\theta$ )} {Thickness or Diameter [cm]}
```

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```
// Total Mass Attenuation Coefficient [cm2/g]
```

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The screenshot shows the 'BL04B2\_Ana' software window. The interface is organized into several panels with different background patterns:

- Top Left (Green background):** Contains a red 'BL04B2 Analysis' header and a 'Load Data' button. Below are input fields for 'Density of sample',  $\mu/\rho$ , 'Polarization Factor', 'B.G. a', 'B.G. b',  $2\theta_0$ , and 'Dead time'. A 'Pol. corr.' checkbox is checked. A button labeled 'Abs, Pol & BKG Correction' is present.
- Middle Left (Green background):** Features a 'Smoothing' checkbox and a 'Combine Data' button.
- Bottom Left (Green background):** Includes 'Auto norm.' (checked), 'Calculation of S(Q)', 'Norm. factor', 'Recoil factor', 'Dump coeff.', and 'fitting to Scattering factor with Q range between 15 and 20'. A 'Subtract Compton scattering' checkbox is checked.
- Bottom Left (Light Green background):** Contains 'Set S(Q=0)', 'Q-range from 0 to 20', and 'with delta-Q 0.05'. A button labeled 'Interpolation of S(Q)' is at the bottom.
- Top Right (Light Green background):** Titled 'Calculation of G(r)', it has a 'Window Function' dropdown set to 'not use', 'r range from 0.01 to 100', and 'with delta-r 0.01'.
- Middle Right (Light Blue background):** Titled 'Set experimental density', it has 'Manual adjust' selected, an 'Auto adjust' option, a 'Show on graph' button, and a 'Number density' input field.
- Middle Right (Light Blue background):** A button labeled 'Transform G(r) into' followed by a dropdown set to 'T(r)'. Below it is a 'Coord. number fit' section with a 'with baseline fit' checkbox.
- Bottom Right (Light Blue background):** A 'Set g(r) to zero' button with a checked 'First range' checkbox, and 'from 0 to 0' input fields. Below is a 'Transform G(r) to S(Q)' button.
- Bottom Right (Light Blue background):** A 'Set this S(Q) as initial' button, followed by a '100% backtransformed S(Q)' dropdown.
- Bottom Right (Light Blue background):** A 'Save dataset' button, followed by a dropdown set to 'Original Q, S(Q) -> RMC'.

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